

Membrane Filters

BF	Biofil PES	F	Fluorofil
C	Chemifil	H	Hydrofil

Absolute Pore Rating

04	0.04µm	20	0.2µm	65	0.65µm
10	0.1µm	45	0.45µm	120	1.2µm

Version

B	Flushed, Economy	L	Economy
D	Rinsed, Economy	R	Rinsed
F	Flushed	S	Standard Hard Cage

Length (Nominal)

1	250mm (10")	4	1020mm (40")
2	510mm (20")	5	125mm (5")
3	760mm (30")		

Adapters - for a comprehensive guide please see publication 2536

A	P Code 3	J	P MYS	R	DOE "Knife Edge"
B	P Code 7	K	P Code 2	S	S Code 28 (No Lugs)
C	Code 8	L	DH Code A (No Lugs)	T	DH Code A (No Lugs)
F	N SOE	M	DOE	V	DH Code C
G	G DOE (Short Length)	N	P Code 3		
H	G SOE	P	P Code 18 (Retro Fit)		

Seals

A	Ethylene Polypropylene	D	Nitrile	G	PTFE Encapsulated Silicone
B	Silicone	E	PTFE Encapsulated Viton		
C	Viton	F	Polypropylene Felt		

Additional Options

A	N+U	P	Pharmaceutical Grade
N	Non Steamable (No Insert)	U	Unbranded

BF20 S2B BP

Example Part Number for a Membrane Filter

Biofil, 0.2µm, Standard Hard Cage, 250mm (10") Long, P Code 7, Silicone Seals, Pharmaceutical Grade

Note: DOE - Double Open Ended, SOE - Single Open Ended



A Division of the Porvair Filtration Group

Fareham Industrial Park, Fareham, Hampshire PO16 8XG UK

Tel: +44(0)1329 285616

Fax: +44(0)1329 822442

Email: info@porvairfiltration.com

Web: www.porvairfiltration.com



FM 00374