

Hydrofil™

Nylon 6.6 Membrane Cartridge Filters



Microbially rated cartridge filters featuring the latest developments in membrane technology, Hydrofil™ cartridges, are based on a naturally hydrophilic nylon membrane.

Hydrofil™ cartridges exploit the narrow pore size distribution and high void volume of the media to provide a choice of cartridges capable of meeting the requirements of most applications. Careful media selection ensures that Hydrofil™ cartridges are very suited to critical particle control down to 0.01 micron ratings. These cartridges offer high flux rates and low differential pressures, a feature common to nylon membranes.

Ordering Information

Product Code: 1 2 3 4 5 6 7

1: Membrane		2: Pore rating		3: Version		4: Length (Nominal)		5: End Fitting		6: Seals		7: Additional	
HT	Hydrofil™	10	0.1µm	R	Rinsed	1	10" (254mm)	A	Code 3	A	Ethylene Propylene	A	N+U
		20	0.2µm	S	Standard	2	20" (508mm)	B	Code 7	B	Silicone	N	Non-steamable (no insert)
		45	0.45µm			3	30" (762mm)	C	Code 8	C	Viton®	P	Pharma Grade
						4	40" (1016mm)	F	N SOE	D	Nitrile	U	Unbranded
						5	5" (125mm)	G	G DOE (short)	E	FEP Encap. Viton®		
								H	G SOE	G	FEP Encap. Silicone		
								J	216 (218), fin	J	DOE PTFE		
								K	Code 2				
								L	223, fin (no lugs)				
								M	DOE				
								S	Code 28, fin (3 lugs)				
								T	223, flat (no lugs)				
								U	224, fin				
								V	226, fin				
								W	F20 +Code 7 (SS Core)				
								X	F20 +Code 2 (SS Core)				
								Y	BS832, flat				
								Z	F20 +Code Y (SS Core)				

Hydrofil™ cartridges benefit from high protein binding characteristics of nylon membranes. They are highly resistant to integrity failure caused by steam sterilisation and have excellent chemical compatibility characteristics. Hydrofil™ cartridges provide a combination of features and benefits previously unavailable from cartridges based on PVDF, mixed esters of cellulose or polysulphone membranes.

Typical Applications

- Biopharmaceuticals: Bioburden reduction and clarification
- Electronics and semiconductors
- Fine chemicals
- Beverages
- Pure water supply (18MΩ.cm)

Features and Benefits

- Guaranteed microbial ratings
- Excellent chemical compatibility
- Cartridge integrity and low TOC levels
- Suitable for steam sterilising
- Full traceability
- Controlled manufacturing environment

Specifications

Materials of Manufacture

Filter membrane:	Nylon 6,6
Membrane support:	Polypropylene
Irrigation mesh (support):	Polypropylene
Drainage layer:	Polypropylene
Inner core:	Polypropylene
Outer support:	Polypropylene
End fittings:	Polypropylene
Support ring:	Stainless steel

Cartridge Dimensions (Nominal)

Effective Filtration Area:	0.63m ² (6.8ft ²) per 10" module	
Diameter:	70mm (2.8")	
Length:	1 module:	254mm (10")
	2 modules:	508mm (20")
	3 modules:	762mm (30")
	4 modules:	1016mm (40")

Other size formats (including juniors) are available upon request.

Cartridge Treatment

Standard:	Cleaned and flushed with pyrogen-free water
Rinsed:	Ultra-clean, pulse flushed to give a system resistivity of 18MΩ.cm

Gaskets and O-Rings

FDA approved Ethylene Propylene, FEP encapsulated, Silicone, Viton® or Nitrile

Maximum Differential Pressure

Normal flow direction at:	
20°C (68°F):	6.0bar (87psi)
80°C (176°F):	4.0bar (58psi)
100°C (212°F):	3.0bar (44psi)
120°C (248°F):	2.0bar (29psi)
Reverse flow direction at:	
20°C (68°F):	2.1bar (30psi)
80°C (176°F):	1.0bar (15psi)
100°C (212°F):	0.5bar (7psi)

Operating Temperature

Maximum continuous: 60°C (140°F)

Sterilisation

In situ steam 40 x 25 min cycles at 121°C (250°F).

Extractables


Minimum total extractables. Please refer to the Hydrofil™ Validation Guide.

Integrity Testing

Each Hydrofil™ module of every cartridge is individually integrity tested using the Diffusive Flow Test, which correlates to the HIMA and ASTM F838-05 bacterial challenge tests. Non-destructive integrity tests, such as Pressure Hold, Diffusive Flow and Bubble Point, can be performed by customers. Please contact us for procedural details.

Clean Water Flow Rates

- Typical clean water flow rate:
A 254mm (10") Hydrofil™ single cartridge exhibits the flow- Δ P characteristics indicated below, for solutions with a viscosity of 1 centipoise.
- Other solutions:
For solutions with a viscosity other than 1 centipoise, multiply the indicated differential pressure by the viscosity in centipoise.



PFG723/Rev12:Oct22